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How to turn on Olympus FV1000 system

1. Turn concentrated power supply ① ON.
2. Turn the LD559nm laser(②)ON
Main (②-1)ON (PL light **Green ON**, TEMP light **Green blinking**).
↓
Wait for TEMP light **Green ON**
↓
Turn the Key(②-2)to the right. LASER light **Red flashes, please wait for Red light ON.**
Please start acquisition after **Red light ON.**
3. Turn on Multi Ar-laser (③)
Main (③-1)ON, then turn the Key(③-2)to the right.
4. Turn the mercury lamp burner (④) ON.
5. Turn the computer ON(⑤), then log on Windows. User name: **fluoview**, password: **fluoview**
6. Double click the icon FV10-ASW to log on to ASW. Username and password: Administrator

How to close Olympus FV1000 system

1. Exit the FV10-ASW software by selecting File/Exit.
2. Exit the Windows
3. Turn the LD559 laser OFF. Turn the Key(②-2) to the OFF position, then MAIN (②-1) OFF.
4. Turn the Multi Ar-laser OFF. Turn the Key(③-2) to the OFF position, then MAIN (③-1) OFF.
5. Turn the mercury lamp burner (④) OFF.
6. Turn concentrated power supply (①) OFF.

※ When finished viewing with 60x oil immersion objective lens, clean the oil immersion objective lens. Use swab moistened with lens cleaning solution. To remove any remaining immersion oil form objective lens, wipe the surface of lens more than once.

Caution

Please restart the system after 30 minutes or more.

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