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



Spatial Omics is Hot, Rapidly emerged field

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In his view, researchers will want single-cell resolution, which scRNA-seq provides but tissue imaging, as of yet, does not. Resolution matters because in a region of interest tissues can have a mix of cell types but "a completely wrong picture" can emerge

CODEX is the leading Spatial Phenotyping solution for Global Cell Atlas Efforts

Global program to create cell atlas


Top tier institute participation


Pacific Northwest
NATIONAL LABORATORY

UC San Diego

Nearly 100 datasets contributed to public domain
Over 50 publications

co-detection by indexing (CODEX)の原理とワークフロー

- Unique Barcode on every antibody
- Complementary oligo attached to each dye enables highly specific detection of Reporters to Barcodes
- Spectrally distinct dyes avoid signal overlap

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
Filter	Excitation (nm)	Emission (nm)	凍結切片	FFPE
DAPI Filter – Nuclear Stain	358	461	✓	✓
AF488/FITC/GFP	492	517	✓	✗
Atto550/TRITC	560	575	✓	✓
Cy5	648	668	✓	✓
AF750 (Cy7)	753	775	✗	✓

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

More flexible in cell
segmentation options →

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Public open source

bioRxiv

Evaluation of cell segmentation methods without reference segmentations

Haoran Chen, Robert F. Murphy

doi: <https://doi.org/10.1101/2021.09.17.460800>

Table 3 Segmentation methods evaluated

Method	Inputs					Output
	Cytoplasm	Cell Membrane	Nucleus	Dimension	Requires scale	Nuclear mask
DeepCell 0.6.0		X	X	2D	No	Yes
DeepCell 0.6.0	X		X	2D	No	Yes
DeepCell 0.9.0 (new)		X	X	2D	No	Yes
DeepCell 0.9.0 (new)	X		X	2D	No	Yes
Cellpose 0.0.3.1	X		X	2D/3D	No	Yes
Cellpose 0.6.1 (new)	X		X	2D/3D	No	Yes
CellSegm		X		2D/3D	Yes	No
CellIX		X		2D	Yes	No
CellProfiler		X	X	2D	Yes	No
AICS (classic)		X	X	2D/3D	Yes	No
Voronoi			X	2D	No	Yes

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


Phenotype Assignment	Population	Positive Markers (All CD45+)
Stromal Cells	3414	CD90+, CD9+, HLA-DR+, Cytokeratin+
Stromal Cells in connective regions	807	CD31+, CD90+, Collagen IV+
T _H Cells	328	CD3+, CD4+, CD7+, CD278+, CD279+
Activated T _H Cells	1459	CD3+, CD4+, CD40+, CD278+, CD279+
Activated Helper T cells	11866	CD3+, CD4+, CD7(+)
T Cells near connective regions	222	CD3+, CD29+, CD90(+), Collagen IV(+)
Proliferating Th-1 cells	2411	CD4+, Ki67+
Dendritic Cells	135	CD4+, HLA-DR+
Dendritic Cells	366	CD4+, CD11c+
Dendritic Cells	429	CD11c+, HLA-DR+
NK cells	126	CD3-, CD7+
Plasma Cells	2180	CD31+, CD38+, CD40+
Lymphatic Vasculature Cells	1336	CD9+, Podoplanin+
Mature, proliferating B cells near follicular matrix	960	CD19+, CD22+, CD21+, CD90+, Ki67+
Mature B cells	32840	CD19+, CD22+

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
- Plasma cells
- Plasma cells near epithelial cells
- DCs near B cells
- NKs
- T cells interacting with DCs
- Vascular endothelial cells (vWF/CD146)
- Combined T cell population
- Epithelial cells (CD63/Pan-cytokeratin)
- Vascular endothelial cells (CD146/CD36)
- Vascular endothelial cells (CD146/vWF/CD36)
- B cells
- DCs interacting with T cells
- Epithelial cells (CD9/vimentin/Pan-cytokeratin)
- Epithelial cells (CD63/CD9/Pan-cytokeratin)
- T cells interacting with DCs
- Helper T cells
- Combined T cell population
- T cells interacting with DCs
- TH1 (CD279+)
- Lymphatic vasculature
- NKs
- Plasma cells
- Plasma cells near epithelial cells
- Vascular endothelial cells (vWF/CD146)
- Vascular endothelial cells (CD146/CD36)
- Filtered

Cellular Neighborhood解析 (CytoMAP)



<https://www.akoyabio.com/webinar/codex-with-seurat-a-high-dimensional-approach-to-unlocking-spatial-insights/>

Cellular Neighborhood解析

Cluster spatial windows to find common patterns of proximity, i.e., cell neighborhoods

- ● ● ● ● ● ● ● ● ● Cell 1 neighbors (window 1)
- ● ● ● ● ● ● ● ● ● Cell 2 neighbors (window 2)
- ⋮

Cellular Neighborhood Interaction

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




Cell

Resource

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



Coordinated Cellular Neighborhoods Orchestrate Antitumoral Immunity at the Colorectal Cancer Invasive Front

Cell. 2020 Oct 29;183(3):838.

nature protocols PROTOCOL
<https://doi.org/10.1038/s41594-021-00556-8>

CODEX multiplexed tissue imaging with DNA-conjugated antibodies

Sarah Black^{1,2,6}, Darci Phillips^{1,2,6}, John W. Hickey^{1,2,6}, Julia Kennedy-Darling^{1,4}, Vishal G. Venkataraman^{1,2}, Nikolay Samusik^{1,5}, Yury Goltsev^{1,2}, Christian M. Schürch^{1,2,3,7,8,9} and Garry P. Nolan^{1,2,7,8,9}

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

Model of differences in the IME between CLR and DII patients with respect to CN organization, cellular function, and inter-CN communication

New Biomarkers are Needed to Better Predict Patient Response

ARTICLE

<https://doi.org/10.1038/s41467-021-26974-6>


OPEN

Check for updates

Immune cell topography predicts response to PD-1 blockade in cutaneous T cell lymphoma

Darci Phillips^{1,2,3,9}, Magdalena Matusiak^{3,9}, Belén Rivero Gutierrez³, Salil S. Bhate^{1,3,4}, Graham L. Barlow^{1,3}, Sizun Jiang^{1,3,5}, Janos Demeter¹, Kimberly S. Smythe⁶, Robert H. Pierce⁶, Steven P. Fling⁶, Nirasha Ramchurren⁶, Martin A. Cheever⁶, Yury Goltsev^{1,3}, Robert B. West³, Michael S. Khodadoust^{7,10}, Youn H. Kim^{2,7,10}, Christian M. Schürch^{1,3,8,10} & Garry P. Nolan^{1,3,10}


Conventional Biomarkers Fail to Predict Outcomes in Lymphoma

Context Matters: Spatial Neighborhoods Predict Response



Discovery Uncovers Predictive Potential of Spatial Neighborhoods

CITN-10 Phase II Clinical Trial (Cutaneous T-Cell Lymphoma)

- 56-plex IO Panel
- Multiple spatial neighborhoods differentiate responders vs. non-responders!

Spatial Neighborhood Map Responder

Spatial Neighborhoods in Today's Literature

Rapid Adoption of SNs as Potential Biomarkers Across the Biomedical Literature

Cell Resource

Coordinated Cellular Neighborhoods Orchestrate Antitumoral Immunity at the Colorectal Cancer Invasive Front

Graphical Abstract

Technology
CODEX + Optimization for Formalin Fixed Paraffin Embedded Tissue (FFPE)

Data Collection
Colorectal Cancer (High Risk, Low Risk) → FFPE → Tissue Microarray

Computation
Cell Types, Cell Neighborhoods
① Coupling
② Functional State
③ Communication

Authors
Christian M. Schürch, Salil S. Bhate, Graham L. Barlow, ..., Nikolay Samusik, Yury Goltsev, Garry P. Nolan

Correspondence
christian.m.schuerch@gmail.com (C.M.S.), gnolan@stanford.edu (G.P.N.)

In Brief
A multiplexed tissue imaging and computational analysis framework applied to colorectal cancer allows interrogation of how spatial organization of the immune tumor microenvironment is linked to clinical outcomes.

Highlights

- FFPE-CODEX multiplexed tissue imaging of 56 markers in 140 tissues of 35 CRC patients
- Cellular neighborhoods reveal spatial organization of the tumor microenvironment

nature COMMUNICATIONS

ARTICLE

<https://doi.org/10.1038/s41467-021-29303-7> OPEN

An N-Cadherin 2 expressing epithelial cell subpopulation predicts response to surgery, chemotherapy and immunotherapy in bladder cancer

Kenneth H. Gouin III^{1,2,8}, Nathan Ing^{1,2,8}, Jasmine T. Plummer^{1,2}, Charles J. Rosser^{3,4}, Bassam Ben Cheikh^{1,2}, Catherine Oh^{1,2}, Stephanie S. Chen^{1,2}, Keith Syson Chan^{4,5}, Hideki Furuya^{3,4}, Warren G. Tourtellotte^{1,4,5,6,7}, Simon R. V. Knott^{1,2,4,8} & Dan Theodorescu^{3,4,5,8}

Cell Systems Article

Tissue schematics map the specialization of immune tissue motifs and their appropriation by tumors

Graphical abstract

Tissue Schematics
Cellular Neighborhoods (CNs): regions with characteristic local processes

Spatial Contexts
Fluctuations in cell states indicate local CN interactions

Assembly Rules
Align tissue types, Evaluate mutations (Good survival, Poor survival)

Authors
Salil S. Bhate, Graham L. Barlow, Christian M. Schürch, Garry P. Nolan

Correspondence
gnolan@stanford.edu

In brief
We present Tissue Schematics, a computational and conceptual framework for using high-parameter imaging data to map how tissues are formed by the assembly of cellular neighborhoods into tissue motifs. Applying our framework to human lymphoid tissues and the colorectal cancer immune-tumor microenvironment suggests mechanistic insights into the specific biology of these tissues, as well as general principles of tissue organization.

nature COMMUNICATIONS

ARTICLE

<https://doi.org/10.1038/s41467-021-26974-0> OPEN

Immune cell topography predicts response to PD-1 blockade in cutaneous T cell lymphoma

Darci Phillips^{1,2,3,9}, Magdalena Matusiak^{3,9}, Belén Rivero Gutierrez³, Salil S. Bhate^{1,3,4}, Graham L. Barlow^{1,3}, Sizun Jiang^{1,3,5}, Janos Demeter¹, Kimberly S. Smythe⁶, Robert H. Pierce⁶, Steven P. Fling⁶, Nirasha Ramchurren⁶, Martin A. Cheever⁶, Yury Goltsev^{1,3}, Robert B. West³, Michael S. Khodadoust^{7,10}, Youn H. Kim^{2,7,10}, Christian M. Schürch^{1,3,8,10,8} & Garry P. Nolan^{1,3,10,8}

Blood Cancer Journal

ARTICLE OPEN

Lack of intrafollicular memory CD4+ T cells is predictive of early clinical failure in newly diagnosed follicular lymphoma

Patrizia Mondello^{1,2,11}, Angelo Fama^{1,11}, Melissa C. Larson¹, Andrew L. Feldman⁵, Jose C. Villasboas¹, Zhi-Zhang Yang¹, Ilya Galkin¹, Viktor Svelolikin¹, Ekaterina Postovalova¹, Alexander Bagaev¹, Pavel Ovcharov¹, Arina Varlamova¹, Sarah Huet¹, Bruno Tesson¹, Kaitlyn R. McGrath¹, Susan Slager¹, Brian K. Link¹, Sergei Syrbu¹, Anne J. Novak¹, Thomas M. Habermann¹, Thomas E. Witig¹, Grzegorz S. Nowakowski¹, Gilles Salles^{1,12}, James R. Cothan¹ and Stephen M. Ansell^{1,13}

CSH Cold Spring Harbor Laboratory bioRxiv

High Resolution Single Cell Maps Reveals Distinct Cell Organization and Function Across Different Regions of the Human Intestine

John W. Hickey^{1*}, Winston R. Becker^{2*}, Stephanie A. Nevins^{2*}, Aaron Horning², Almudena Espin Perez³, Roxanne Chiu², Derek C. Chen², Daniel Cotter², Edward D. Esplin², Annika K. Weimer², Chiara Caraccio¹, Vishal Venkataaraman¹, Christian M. Schürch^{1,7}, Sarah Black¹, Maria Brbic⁵, Kaidi Cao⁵, Jure Leskovec⁵, Zhengyan Zhang⁵, Shin Lin⁴, Teri Longacre¹, Sylvia K. Plevitis³, Yiing Lin⁵, Garry P. Nolan¹*, William J. Greenleaf²*, Michael Snyder²*

Footnotes:
¹ Department of Pathology, Stanford School of Medicine, Stanford, CA, 94305
² Department of Genetics, Stanford School of Medicine, Stanford, CA, 94305
³ Department of Biomedical Data Science, Stanford School of Medicine, Stanford, CA, 94305
⁴ Department of Medicine, University of Washington, Seattle, WA, 98195
⁵ Department of Surgery, Washington University, St Louis, MO 63110
⁶ Department of Computer Science, Stanford University, Stanford, CA 94305, USA
⁷ Department of Pathology and Neuro pathology, University Hospital and Comprehensive Cancer Center Tübingen, Tübingen, Germany

*These authors contributed equally
Correspondence should be addressed to: gnolan@stanford.edu, wjg@stanford.edu, mpsnyder@stanford.edu

Research Square

Preprints are preliminary reports that have not undergone peer review. They should not be considered conclusive, used to inform clinical practice, or referenced by the media as validated information.

A tissue atlas of ulcerative colitis to guide TNF inhibitor therapy


Aaron Mayer, Derek Holman, Utkarsh Tandon, Salil Bhate, Sunil Bodapati, Graham Barlow, Jeff Chang, Sarah Black, Erica Crenshaw, Brigid Boland, Sarah Streett, William Sandborn, John Chang, Sanjiv Gambhir, Garry Nolan, Christian Schürch, Stephan Rogalla

Abstract
Ulcerative colitis is a chronic-relapsing inflammatory disease of the large intestine with a complex, multifactorial pathogenesis. TNF inhibitors are widely used to suppress immune-mediated tissue damage in ulcerative colitis patients; however, therapy failures are common. Predicting TNF inhibitor response requires an understanding of the architectural features that underlie mucosal inflammation and those responsible for resistance. Here, we used highly multiplexed immunofluorescence to uncover the spatially resolved tissue architectures underlying disease progression and treatment response in 42 tissue regions from 34 individuals. We created a tissue atlas and performed spatial analysis to identify cell-cell contacts and cellular neighborhoods. We observed that cellular functional states depend on cellular neighborhood and that a subset of inflammatory cell types and cellular neighborhoods in ulcerative colitis patients persisted even during treatment with TNF inhibitor, indicating resistant niches. A computer vision model, with no *a priori* assumptions regarding cellular architectural features, was able to predict TNF inhibitor resistance. This spatial model significantly outperformed classification models based on single-cell data alone. Our results demonstrate the value of a spatial tissue atlas as a precision medicine tool to guide treatment of patients suffering from autoimmune diseases.

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

CODEX, Together with Sn-Sequencing, Identifies Bio-marker for Cancer

An N-Cadherin 2 expressing epithelial cell subpopulation predicts response to surgery, chemotherapy and immunotherapy in bladder cancer

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CODEX : True Single Cell Level Protein Expression With Spatial Information


JCI insight

Integration of spatial and single-cell transcriptomics localizes epithelial cell-immune cross-talk in kidney injury

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

Multiplexed imaging of proteins in toto with **CODEX** validates the localization of immune cell clusters inferred by Spatial Transcriptomics

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Merger of CITE-Seq and PhenoCycler data

CITE-Seq produces multiplexed PROTEIN and mRNA expression data

CITE-Seq workflow

PhenoCycler (CODEX) results are multiplexed PROTEIN expression data

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
Resource

Deep Profiling of Mouse Splenic Architecture with CODEX Multiplexed Imaging

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
ScienceAdvances

Single-cell transcriptomic analysis of mIHC images via antigen mapping

Science Advances 05 Mar 2021:
Vol. 7, no. 10, eabc5464
DOI: 10.1126/sciadv.abc5464

Kiya W. Govek*, Emma C. Troisi*, Zhen Miao, Rachael G. Aubin, Steven Woodhouse, Pablo G. Camara†

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Human Biomolecular Atlas Program (HuBMAP) 関連

Highly Multiplexed Immunofluorescence of the Human Kidney using CODEX

Kidney International

Volume 101, Issue 1, January 2022, Pages 137-143


Technical Notes

Highly multiplexed immunofluorescence of the human kidney using co-detection by indexing

Elizabeth K. Neumann^{1,2}, Nathan Heath Patterson^{1,2}, Emilio S. Rivera^{1,2}, Jamie L. Allen^{1,2}, Maya Brewer³, Mark P. deCaestecker³, Richard M. Caprioli^{1,2,4}, Agnes B. Fogo^{3,5,6,7}, Jeffrey M. Spraggins^{2,4,8} ✉

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CellPress
OPEN ACCESS

Protocol

Protocol for multimodal analysis of human kidney tissue by imaging mass spectrometry and CODEX multiplexed immunofluorescence

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



TABLE 1	SPECIMEN	
Tissue Type	FFPE Oropharyngeal Squamous Cell Carcinoma	FFPE Human Tonsil
Markers	103 Antibodies	103 Antibodies
Tissue imaging area	1.2 cm x 1.4 cm	1.3 cm x 1.8 cm
Fluidics time	38 cycles 44 hours	38 cycles 50 hours
Data size	24 gigabytes	46 gigabytes
Number of cells	832,823 cells	2,253,539 cells

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RNAScope HiPlex v2

Automated HiPlex v2 Workflow on PhenoCycler-Fusion System

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3. Spatial Transcriptomics Workflow is Compatible with a Broad Set of Formalin-fixed, Paraffin-embedded (FFPE) Tissues

3.1 104-Plex RNA Workflow on Human Tissue Microarray

RNA detection across a human tissue microarray on the PhenoCycler-Fusion with a 104-plex RNA panel. Detection of RNA targets produces distinct labeling patterns in different tissue types. All cores are 2mm in diameter and all samples are human FFPE.

3.2 102-Plex RNA Workflow on Human HNC Squamous Cell Carcinoma Sample

RNA spatial phenotyping of a human head and neck cancer tissue on the PhenoCycler-Fusion with a 102-plex RNA panel. The images show a whole slide view (left) and zoomed-in ROIs to illustrate different marker combinations and tissue features.

5. Assay Compatibility Enables Same-Section Multiomic Spatial Phenotyping

Same section, spatial multiomic data

ACKR1 CD21-P CD45-P CD68 ICHAM1 MPO-P PanCK-P

B

- B cells
- Connective Tissue
- Endothelial cells
- Epithelial cells
- Glandular cells
- Macrophages
- Necrotic Tissue
- Other
- Plasma cells
- T cells
- Tumor cells


C

CD68-P ANXA1 IML KRT10 DAPI TP53 Keratin 5-P

D

CD68-P ANXA1 IML KRT10 DAPI TP53 Keratin 5-P

(A) Same-section, multiomic spatial phenotyping of head and neck cancer tissue. The compatibility of our high-plex protein and RNA detection assays enables multiomic analyses on the same tissue section. (B) Spatial projection of key cell types identified through clustering of the multiomic dataset. (C) The same ROI as shown in the RNA-only panel in Section 4.1, but with protein targets highlighted for emphasis. (D) The same ROI as shown in the RNA-only panel, but with additional protein markers, including additional Keratins and CD68 to denote macrophages. This structure shows the attack of a necrotic region of tissue and potential destruction of the Keratin border.

Spatial Multiomics Grant Program
Rapid, Unbiased *In Situ* RNA & Protein Detection
 Application Deadline: September 30, 2022

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